# FX CorAL

## HD Dialyzers Single Use Only – INLINE Steam Sterilized

**Technical Data:** These data represent typical *in vitro* performance measured according to ISO 8637-1. Actual *in vivo* performance may differ due to the patient's blood composition and clinical settings. Performance values should be regarded as approximate. Results may vary due to differences in measurement method as well as batch variation.

Model		FX CorAL 60	FX CorAL 80	FX CorAL 100	FX CorAL 120
Catalogue number		F00009216	F00009217	F00009218	F00009219
Membrane Surface Area	(m <sup>2</sup> )	1.4	1.8	2.2	2.5
Blood Ultrafiltration Coefficient <sup>1</sup>	(mL/hr/mmHg)	61	75	84	94
Priming Volume Blood	(mL)	80	102	124	142
K <sub>0</sub> A Urea <sup>3</sup>	(mL/min)	1145	1321	1487	1614
Maximum pressure	(bar)	3			
Pressure drop	(mmHg)				
Across blood compartment	(Q <sub>B</sub> = 300 mL/min)	85	66	57	48
Across dialysate side	(Q <sub>D</sub> = 500 mL/min)	41	34	27	23
Maximum TMP	(mmHg)	600			
Minimum Blood Flow	(mL/min)	150	200	250	300
Maximum Blood Flow	(mL/min)	400	500	600	600
Minimum Dialysate Flow	(mL/min)	300	300	500	500
Maximum Dialysate Flow	(mL/min)	1000	1000	1000	1000
Hemodialysis (HD) Clearance <sup>2</sup> (mL/min) Q <sub>D</sub> = 500 mL/min, Q <sub>F</sub> = 0 mL/min	Q <sub>B</sub> (mL/min)	200 300 400	200 300 400	300 400 500	300 400 500
	Urea <sup>3</sup>	195 270 317	195 277 328	282 339 372	285 345 381
	Creatinine	190 246 281	195 256 293	264 305 332	272 315 343
	Phosphate	188 248 287	191 255 297	266 313 343	268 320 351
	Vitamin B <sub>12</sub>	143 167 183	155 184 201	197 216 230	209 233 248
	Inulin	88 99 102	101 110 118	123 133 136	129 139 146
	Cytochrome C <sup>4</sup>	81 84 86	95 99 103	117 118 121	126 128 129
Sieving Coefficient <sup>1,5</sup>	Albumin	< 0.001			
	Myoglobin	0.6			
	β <sub>2</sub> -Microglobulin	1			
	Inulin	1			


For full range of clearance and sieving results see:  
<https://fmcna.com/products/disposables/dialyzers/fx-coral>

### Dialyzer Technical Information

Membrane Material	Helixone <sup>®</sup> hydro
Fiber Inner Diameter	185 μm
Membrane Wall Thickness	35 μm
Housing Material	Polypropylene
Potting Compound	Polyurethane
O-ring	Silicone
Protective Caps	Styrene-ethylene-butylene-styrene (SEBS)
Recommended Blood Connections	ISO 8637-1:2017
Recommended Dialysis Fluid Connections	ISO 8637-1:2017
Sterilization Method	INLINE Steam

### Footnotes:

- <sup>1</sup> Human blood, Hct 32 %, protein 60 g/L, 37 °C.
- <sup>2</sup> Aqueous solution of respective test substance, 37 °C.
- <sup>3</sup> Sodium is used as a surrogate for urea.
- <sup>4</sup> Cytochrome C (12.3 kDa) is used as a marker for middle-molecule clearance.
- <sup>5</sup> Average of sieving coefficients measured with human blood at Q<sub>B,max</sub> and Q<sub>B,min</sub>. In both cases, the filtrate flow is 20 % of Q<sub>B</sub>.

**Figure 1:** Position of Dialyzer during Priming and Treatment



**FRESENIUS  
MEDICAL CARE**

Manufactured for:  
**Fresenius Medical Care Renal Therapies Group LLC**  
Waltham, MA 02451  
☎ 1-800-323-5188

# FX CorAL HD Dialyzers

## GENERAL INFORMATION

Refer to product or carton label for:

	Single use only		Blood pathway sterile. Steam sterilized
	Expiry date		Fluid path is non-pyrogenic
	Batch code		Refer to instructions for use
	Date of manufacture		Caution
	Catalogue number		Temperature limit
	Units		Caution: Federal law restricts this device to sale by or on the order of a physician.

**Indications for Use:** FX CorAL HD dialyzers are intended for hemodialysis (HD), hemodiafiltration (HDF), hemofiltration (HF) and isolated ultrafiltration in patients with acute kidney injury or chronic kidney disease when conservative therapy is judged to be inadequate.

**Contraindications:** FX CorAL HD dialyzers are not intended for use in patients that weigh less than 30 kg.

Patients with known hypersensitivity to any of the materials used (see Dialyzer Technical Information) must not be treated with this product.

There are no absolute contraindications to hemodiafiltration; however, hemodiafiltration should be avoided or applied with extreme caution in patients with the following conditions or circumstances:

1. Unable to receive or contraindications to adequate anticoagulation.
2. Hypercoagulation syndromes.
3. Hyperviscosity syndromes.
4. Pre-treatment hematocrit  $\geq 40\%$
5. Single needle blood access.
6. Blood access with obvious recirculation.

**Dialysate:** The dialysate must meet current water quality standards for dialysis (AAMI RD52 and ISO 11663) to mitigate risk of infection due to back-filtration.

**Side effects:** Adverse reactions such as hypertension, hypotension, headache, and nausea may be associated with hypovolemia or hypervolemia and can usually be avoided with careful management of patient fluid, electrolyte balance, blood flow rate, and ultrafiltration rate. During treatment adequate time for administration of accompanying medication should be considered.

In rare cases, thrombocytopenia or hypersensitivity reactions including anaphylactic or anaphylactoid reactions to the dialyzer, or other elements in the extracorporeal circuit, may occur during treatment. Hypersensitivity reactions may cause mild to severe signs and symptoms, including: itching, flushing, hives, swelling, fever, leukopenia, hypotension, hypertension, shortness of breath with wheezing, arrhythmias, and/or respiratory arrest.

Patients with a history of hypersensitivity reactions or patients who have a history of being highly sensitive and allergic to a variety of substances should be carefully monitored during treatment. With severe hypersensitivity reactions, dialysis must be discontinued and aggressive first line therapy for hypersensitivity reactions must be initiated. Blood from the extracorporeal system should not be returned to the patient in cases of hypersensitivity reactions as determined by the physician.

The dialyzer is steam sterilized and thus contains no sterilization residues. For patients with a known ethylene oxide hypersensitivity all components of the extracorporeal circuit should be ethylene oxide-free.

Hemolysis may also occur as a result of dialysis.

**Anticoagulation:** We recommend systemic anticoagulation of the patient. Systemic anticoagulation is the administration of the prescribed loading dose of anticoagulant into the vascular access of the patient; then waiting a period of 3 to 5 minutes prior to initiating the treatment. During dialysis, the dose of anticoagulant and method of administration is the decision of the physician.

**Sterile/Non-pyrogenic:** The blood pathway is sterile and non-pyrogenic if blood port caps are in place and undamaged. Do not use if the dialyzer is damaged in any way. Use aseptic technique for all blood side connections. Structural integrity of the dialyzer is warranted when prepared as directed. Sterilized using the INLINE steam method of sterilization.

**Dialyzer reuse:** FX CorAL HD dialyzers are intended for single use only. Re-use of the device may create a potential risk to the patient, including contamination and/or impairment of the device function.

## PRECAUTIONS

Dialyzers may leak resulting in patient blood loss or contamination with dialysate. In the event of a blood leak during dialysis, the health care provider should respond according to the facility's established protocol.

Only sterile, nonpyrogenic fluids labeled for direct infusion should be used as replacement fluids during hemodiafiltration and hemofiltration therapies, otherwise serious or adverse events may occur.

Strict patient monitoring is recommended during treatments as there may be a potential for hypotensive episodes.

## CAUTIONS

Federal (US) law restricts this device to sale by or on order of a physician.

The operator should strictly adhere to the manufacturer's recommended procedures, warnings, and cautions as listed in these instructions for use. Appropriate adjustments to drug doses should be made to account for their dialytic removal.

Selection of the appropriate dialyzer should take into consideration underlying patient health conditions. Caution for use of the dialyzer should be considered in certain patients, including patients with low weight or blood volume, who are elderly, a history of hypotension, who are pregnant, with severe heart failure and, on life support.

Provide increased caution if using the device in patients recently initiating haemodialysis.

## WARNINGS

Due to the high water flux capability of high permeability membranes with an ultrafiltration coefficient  $> 8$  mL/hr/mmHg, it is necessary to use such dialyzers only in conjunction with dialysis machines that are equipped with precise ultrafiltration control. In any case, the safety instructions for the hemodialysis machine must be followed. Only use the device in dialysis machines that have intact blood leak detectors and air detectors.

Dialyzer performance may change with the duration of observation.

The user is cautioned to regularly monitor the patient's chemistry values using quantitative measurements and analysis to ensure that the prescribed therapy is delivered. The clinical parameters monitored should at least include urea, hemoglobin and serum albumin.

During postdilution-HDF and postdilution-HF with manual setting of the substitution rate, the maximum filtration rate (= substitution rate + UF rate) must not exceed 25 % of the effective blood flow rate. This is to keep the blood water content within a certain noncritical range and to avoid the risk of hemoconcentration and coagulation in the extracorporeal circuit. Hemoconcentration and coagulation would be reflected by increased TMP values.

Dialyzers must not be used after expiry date (see label).

Blood flow rates below the recommended blood flow range can be set. In this case the coagulation should be monitored by a standard clotting test and the decrease of the clearance values might be steeper than expected.

Blood flow rates above the recommended blood flow range can be set. This might result in an increase of the sieving coefficient of albumin when performing HDF or HF. Further this might result in a lower increase of middle molecule (e.g.  $\beta_2$ -m, myoglobin) clearance than expected.

Do not use the device if package is damaged or if protective end caps are not in place.

Air entering the extracorporeal circuit during dialysis treatment phase can result in air embolism, serious injury or death. Check the security of all extracorporeal connections prior to the initiation of dialysis and periodically at regular intervals throughout the treatment. In the event that air enters the extracorporeal circuit and cannot be removed safely, discontinue the treatment. In such a circumstance, do not return the blood to the patient.

## PREPARATION FOR DIALYSIS

Note: Do not remove blood port caps until ready to attach arterial and venous bloodlines to the dialyzer

- Place the dialyzer in the dialyzer holder in the vertical position.
- Connect the arterial and venous blood lines to the lower and upper dialyzer ports respectively. Avoid contamination of the fluid path of the connection.
- Refer to Figure 1: Position of Dialyzer during Priming and Treatment.
- Attach the dialysate tubing connectors with the inlet to the upper and the outlet to the lower dialysate port.
- Prime the dialyzer according to the instructions provided with the machine (there is no need to rotate the dialyzer). Recommended priming volume: 500 mL. Ensure complete de-aeration of the blood compartment.
- Recirculation: If applicable, refer to the machine operator's manual for instructions on recirculating the bloodline.

## INITIATION OF DIALYSIS

- Ensure that the blood tubing system and dialysate couplings are firmly connected to the dialyzer before the patient is connected.
- Connect the patient to the blood line according to the instructions provided with the machine.

## DURING THE DIALYSIS TREATMENT

- If a blood leak should occur during the treatment, the operator should follow established facility procedures for a dialyzer blood leak.

## COMPLETION OF DIALYSIS

- Perform reinfusion of the blood at the end of the treatment according to the instructions provided with the machine.
- Discard the extracorporeal circuit in an appropriate biohazard waste receptacle. References: 29CFR, 1910.145, 1910.1030 (Code of Federal Regulations) and appropriate state and local codes.

**FRESENIUS  
MEDICAL CARE**

**Manufactured for:**  
Fresenius Medical Care Renal Therapies Group LLC  
Waltham, MA 02451  
☎ 1-800-323-5188

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