



Fresenius Medical Care SERVICE BULLETIN

Copyright © 2007 Fresenius USA, Inc

Equipment: All Hemodialysis Machines
Bulletin: 07-FRS-001 Rev. A
Subject: New Tubing Packaging

Originated By (Dept):	Reviewed By:	Approved By:	Effective Date:
J. Fox (TSD)	Brooks E. Rogers	ECO# 23588	9-27-07

Distribution: Dir TS (B. Rogers), TS Mgr (K. Cogswell), Dir QS (M. Owen), QE (R. Fox), Dir Eng (D. Gavin), DTS (M. Crnkovich), Tng (J. Barr)

1.0 PURPOSE

To inform all Field Service and Technical Support personnel of a change in how some of our tubing will be sold as a spare part.

2.0 SCOPE

As of October 1, 2007, Fresenius will be offering certain types of tubing in pre-packaged lengths. These tubing types will only be available in pre-packaged 10 foot lengths or as complete rolls.

The following table indicates which types of tubing are included in the change as well as the new part numbers.

Description	Package Part Number (1 each = 10 feet)	Roll Part Number
PVC Clear Vinyl w/ nylon braid, 10mm ID x 16mm OD	545210-10 (10 foot package)	545210 (60 meters = 196 feet)
Silicone, Clear, 3mm ID x 7mm OD	545325-10 (10 foot package)	545325 (50 meters = 165 feet)
Silicone, White, w/ Nylon reinforcing, 3mm ID x 6mm OD	545356-10 (10 foot package)	545356 (50 meters = 165 feet)
Silastic, Clear, 5mm ID x 11mm OD	545361-10 (10 foot package)	545361 (50 meters = 165 feet)
Silastic, White, w/ Nylon reinforcing, 6mm ID x 12mm OD	545372-10 (10 foot package)	545372 (50 meters = 165 feet)
Silicone, Clear, w/ Nylon reinforcing, 5mm ID x 11mm OD	N/A	642391 (10 meters = 32.8 feet)

If there are any questions regarding this bulletin, contact Fresenius Technical Support at 800-227-2572 x7003.