



Fresenius Medical Care

SERVICE BULLETIN

Copyright © 2008 Fresenius USA, Inc

Equipment: All Hemodialysis Machines
Bulletin: 08-FRS-001 Rev. A
Subject: Blood Pressure Module Kits (Update)

Originated By (Dept):	Reviewed By:	Approved By:	Effective Date:
J. Fox (TSD)	Kenneth Cogswell	ECO# 25308	June 10, 2008

Distribution: Dir TS (B. Rogers), TS Mgr (K. Cogswell), Dir QS (M. Owen), QE (R. Fox), Eng (D. Gavin), QC Mgr (D. McKenzie), Tng (J. Barr), Prod Mgr (R. Few)

1.0 PURPOSE

To provide all Field Service and Technical Support personnel with a description of parts included in each model Blood Pressure Module Kits used in the 2008H, 2008K and 2008K@HOME hemodialysis machines. This Service Bulletin replaces FSB# 05-FRS-002.

2.0 DESCRIPTION

The following tables indicate the parts included in each kit. Please note that each module type uses a unique cable, fastening hardware and installation instructions. The latest revision of the mounting bracket has all of the appropriate holes to mount all module types.

TM-2910 Kit (P/N 150479)	
370061	Blood Pressure Module TM-2910
150478	Cable, TM-2910
290075	Mounting Bracket
552244	2 Screws, M3x8 SS
330608-02-01	2 Washers, #6 Int Tooth Lock Washer - SS
490070	Quick Reference Card, Blood Pressure Module
507641	Kit Installation Instructions TM-2910

M2600 Kit (P/N 150485)	
370080	Blood Pressure Module M2600
150481	Cable, M2600
290075	Mounting Bracket
330920-02-01	3 Screws, PNP, M3x6 with Lock Washer – Zinc
490070	Quick Reference Card, Blood Pressure Module
507711	Kit Installation Instructions M2600

SunTech Kit (P/N 190474)	
370085	Blood Pressure Module SunTech
190225	Cable, M2000
290075	Mounting Bracket
330611-10	4 Nuts, Nylon, M3 X 0.5
330920-14-01	4 Screws, PH, PHL, with Lock Washer, M3X10,SS
332151-05	Elbow Connector, 5/32X5/32 Barb
330916-10	4 Standoffs Male-Female M3 12MM Long
490070	Quick Reference Card, Blood Pressure Module
507859	Kit Installation Instructions SunTech

If there are any questions regarding this bulletin, contact Fresenius Technical Support at 800-227-2572 x7003.